


# Sample Prep and Experimental Design

Susan Kloet

Leiden Genome Technology Center

11 October 2022




# What is sample prep?



# Know what you are working with

Cells Tested	Species	Cell Source	Total RNA (pg/cell)	Cell Size ( $\mu\text{m}$ )
PBMC	Human	Extracted from blood	~0.75	~5-10
E18 neuron	Mouse	Brain tissue	~2-3	~9
Jurkat	Human	Suspension	5.5	~12
Raji	Human	Suspension	7.3	~12
293T	Human	Adherent	14.2	~18
3T3	Mouse	Adherent	16.1	~18
HCC1954	Human	Adherent	15.7	~18
HCC38	Human	Adherent	21.6	~30

# Quality is critical

<p>Goal is to minimize</p>	Clean	Healthy	Gentle
			
	<ul style="list-style-type: none"><li>• Aggregates/clumps</li><li>• Subcellular debris</li><li>• Free-floating RNA/DNA</li></ul>	<ul style="list-style-type: none"><li>• RNA leakage</li><li>• RNA degradation</li></ul>	<ul style="list-style-type: none"><li>• Physical decomposition</li></ul>

# Handle with care

Minimize handling, both pipetting and centrifugation

Metric	Control	Wide Bore (Rough)	Narrow Bore (Rough)	Vortex 5s
Number of Cells	1,118	846	1,012	983
Reads per Cell	50,000	50,000	50,000	50,000
Valid Barcodes	95.40%	95.50%	95.30%	95.50%
Fraction Reads in Cells	79.40%	72.80%	54.00%	63.10%
Reads Mapped Confidently to Transcriptome	70.50%	71.40%	71.80%	71.00%
Valid UMIs	99.40%	99.40%	99.40%	99.40%
Median Genes per Cell	3,137	3,180	2,833	2,934

# Sample procurement, storage, transport

## CNV

Cultured cells, PBMCs  
\*Fresh and cryopreserved  
\*10x protocols work as-is

Fresh tissue  
\*Minimal optimization

Snap-frozen tissue  
\*Minimal optimization

Fixation?  
\*Not tested but likely  
(ethanol or methanol)

## ATAC

Cultured cells, PBMC  
\*Fresh and cryopreserved  
\*10x protocols work as-is

Fresh tissue  
\*Optimization

Snap-frozen tissue  
\*Optimization

Fixation?  
\*Not tested

## GEX


Cultured cells, PBMC  
\*Fresh and cryopreserved  
\*10x protocols work as-is

Fresh tissue  
\*Significant optimization

Snap-frozen tissue  
\*Nuclei (no cell surface proteins)  
\*Technically challenging  
\*Significant optimization

Fixation?  
\*Methanol (but we prefer  
cryopreservation)

# How to choose a dissociation method



# Separation

## **Separate intact cells and nuclei from**

- Aggregates/clumps
- Debris
- Free-floating mRNA
- Dead cells
- Enrichment/depletion

## **Challenges with separation**

- Samples are fragile
- Physical stress
- Buffers
- Time
- Yield

Want the minimum handling necessary. Maintain sample integrity.





# There are many separation methods

Method	Thorough c'fuge (e.g. 3x in PBS + 0.04% BSA)	Gentle c'fuge (e.g. 1x in media)	Magnetic beads	Density Gradient	FACS
10x Protocol Example	PBMC (CG000039)	Cell Prep Guide (CG000053)	Dead cell removal (CG000093)	Nuclei isolation (CG000124)	Customer Developed Protocol (Martelotto)
Can be used with	Abundant	Limited	Abundant	Abundant	Limited
Benefit	Clean	Gentle	Specific, easily accessible, scalable	Removes Debris	Versatile, quick
Challenges	Lossy, can be harsh	Less thorough	Lossy	Lossy, harsh	Expensive, can be harsh




# Washing / straining

- 3 washes in PBS + 0.04% BSA before loading
- Use a relatively large volume for washes relative to cell number
  - Concentrated cells can clump
- Large cell aggregates or debris can increase the risk of clogs
  - => Use a cell strainer, 30 to 40  $\mu\text{m}$ 
    - BUT – this will result in sample loss, volume loss, and possible cell loss

# Counting and viability assessment




# Counting and viability assessment – live/dead



- Alive
- Dead
- Aggregation


# Dead cell and debris removal

## Debris removal protocol overview, steps 6–13




Source:  
Miltenyi


1.  
The ready-to-use  
MicroBeads are added to  
the cell preparation.




2.  
The sample is separated over  
a MACS® Column in the  
magnetic field of a MACS®  
Separator. Cell debris, dead,  
and dying cells are retained  
in the column.



3.  
The viable cell fraction  
is collected in the flow-  
through. Cells can now  
be used directly for  
culture and functional  
experiments.



# Example: tumor processing workflow



# Summary of key sample prep lessons

- Treat cells gently and keep them happy
  - Try gentle(r) lysis conditions
  - Reduce washing steps
  - Use a swinging bucket centrifuge
  - Keep cells in media + FBS instead of PBS
- Work quickly
  - Consider sorting
  - Minimize handling steps
- *Consider the benefits and drawbacks of different techniques*

# Storage after preparation? Older options


- [Cryopreservation](#) (DMSO)
  - No quantifiable difference between fresh and cryopreserved scRNA-seq samples
- [Methanol fixation](#)

**\*\*NOTE:** both approaches have shown a potential bias in cell-type composition, and it is strongly recommended to thoroughly evaluate preservation methods for new cell types that have not been tested.



# Storage after preparation? New options!

- Fixed RNA kit from 10x Genomics





\*4% PFA fixation  
Storage at 4°C short-term or  
-20°C long-term

- Parse Biosciences (SPLiT-seq)



# Storage after preparation? Other options!

- iCell8 and Rhapsody microwells – freeze plate after single cell sorting
- The HIVE from Honeycomb – freeze HIVE after collection



- Low(er)-throughput option – sort into plates using FACS
  - Great for the DIY workflows like SMART-seq, CEL-seq2, VASA-seq, etc.





# Nuclei vs whole cell

What if cells are not an option?

- Neurons
- Snap-frozen
- FFPE


snRNA-seq is nearly identical to scRNA-seq

- Fewer unique transcripts
- More variability in per-nucleus content



How do I design my single-cell experiment?


# Cell hashing / multiplexing option 1



Pooled 8 PBMC samples in 1 reaction, “super-loading” 20k cells per reaction (5k multiplet)

## Demuxlet

Demultiplex samples using natural genetic variation



# How many cells?

Consider sample heterogeneity and subpopulation frequency

General rule of thumb:

- Heterogenous populations require MORE cells
  - Here, a larger number of cells will provide greater statistical power
  - Heterogeneity can be reduced through experimental design (sorting, transgenics)
- Identification of rare cells requires MORE cells

# How many samples?

## Power calculations + group size estimates


- Is other single-cell data available (flow- or mass-cytometry data)?
- Computational tools available:
  - [How many cells?](#) – Satija lab
  - [SCOPIT](#) – Navin lab
  - [scPower](#) – Theis and Hening labs

Study parameters

Samples	5
Cell type frequency	0.01
Detection power	0.95
Minimal number of cells	<input type="range" value="10 22"/>

1 11 21 31 41 51 61 71 81 91 100

Required cells per person to detect rare cell types with a certain power



The figure shows the required number of cells per individual (y-axis, log scale) to detect the minimal number of cells from a target cell type per individual (x-axis) with a certain probability. The power depends on the total number of individuals and the frequency of the target cell type. Note that the required number of cells per sample only counts correctly measured cells (no doublets etc), so the number is a lower bound for the required cells to be sequenced.

Pilot experiment recommended before large-scale data production




# How deep?

For 3' counting approaches:

<15k reads/cell, doubling reads reduces error by 30-40%

>15k reads/cell, doubling reads or cells has little effect (1-3%)



# How deep?

For 3' counting approaches:


One read per cell per gene

[sceb \(single-cell empirical Bayes\) estimator \(Python package\)](#)

Zhang, Ntranos and Tse Nat Comm 2020


a

Sequencing budget allocation problem



# How deep?

- BUT... not all cells will have the same depth of sequencing



# Conclusions

- Think before you start
  - Which library prep protocol will you choose?
  - How should you store your samples?
  - How should you prepare your tissue/cells/nuclei?
- Maintain sample integrity
- Minimize batch effects
- Make sure your experiment/study has enough power